



Email : sales@laboratoryinstrumentindia.com

Product Name :
Cheese Vat And Cheese Making

Product Code :
LBNY-0005-10700028



Description :

Water heater fitted in console
Temperature and pH can be monitored by microcomputer using an appropriate interface
Ten litre capacity stainless steel jacketed vat
Bench mounting with small footprint
Simple controls
Range of cheese making accessories
Separate console controls vat temperature and paddle speed.

Technical Specification :

Cheese Vat And Cheese Making



Laboratory Instrument India