



Email : sales@laboratoryinstrumentindia.com

Product Name :
HTST UHT Process System

Product Code :
LBNY-0005-10700029



Description :

Full sterile capability options
Touch screen control panel for ease of use
Hygienic fittings as standard
Totally Modular System
Direct and/or indirect heating
Scraped surface, Tubular & Plate heat exchangers
Controllable Pre-heat - option
Built-in CIP facilities
Integral homogeniser option
Standard throughputs from 12 60 l/h
Process temperatures >150C
Electronic flow meter - option
USB data logging - option.

Technical Specification :

HTST UHT Process System

Laboratory Instrument India

