



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
HTST UHT Process System

**Product Code :**  
LBNY-0005-10700029



**Description :**

Full sterile capability options  
Touch screen control panel for ease of use  
Hygienic fittings as standard  
Totally Modular System  
Direct and/or indirect heating  
Scraped surface, Tubular & Plate heat exchangers  
Controllable Pre-heat - option  
Built-in CIP facilities  
Integral homogeniser option  
Standard throughputs from 12 60 l/h  
Process temperatures >150C  
Electronic flow meter - option  
USB data logging - option.

**Technical Specification :**

HTST UHT Process System

---

## Laboratory Instrument India

