




Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Homogenisation Sub Systems

**Product Code :**  
LBNY-0005-10700030



**Description :**

Complete Homogenisation Sub systems  
Comprehensive instrumentation includes product temperature, product line pressure, homogenisation pressure and pump speed.  
Adjustable divert facility for easy integration into continuous flow processes.  
Two Stage Homogenisation with adjustable pressure .  
Twin Pistons or Single Piston  
Pulsation Dampers included  
Configurable for upstream or downstream homogenisation  
Variable throughput  
Can be sterilised and used in an Aseptic system  
Complete with interconnection hoses.

**Technical Specification :**

Homogenisation Sub Systems



**Laboratory Instrument India**

