




Email : sales@laboratoryinstrumentindia.com

Product Name :
Homogenisation Sub Systems

Product Code :
LBNY-0005-10700030



Description :

Complete Homogenisation Sub systems
Comprehensive instrumentation includes product temperature, product line pressure, homogenisation pressure and pump speed.
Adjustable divert facility for easy integration into continuous flow processes.
Two Stage Homogenisation with adjustable pressure .
Twin Pistons or Single Piston
Pulsation Dampers included
Configurable for upstream or downstream homogenisation
Variable throughput
Can be sterilised and used in an Aseptic system
Complete with interconnection hoses.

Technical Specification :

Homogenisation Sub Systems



Laboratory Instrument India

