



Email : sales@laboratoryinstrumentindia.com

Product Name :
Carbonator Filler

Product Code :
LBNY-0005-10700031



Description :

Accurate control of filling and carbonating variables
Throughputs up to 60l/h achievable (2l bottles)
Wide range of cans or bottles can be used from 0.15l to in excess of 2l
Evaluation service for customer cans and bottles
Option of filling and crown seal capping and/or screw top capping at a single station
15l process vessel, positioned at a convenient height
Batch & continuous processing
CO₂ content in excess of 10g/l can be realised
Enhanced colour touch screen control
Carbonates at any temperature between 4°C and 20°C
Premix and postmix capability
Fills non-carbonated drinks
In-bottle pressure measurement enables accurate and reproducible depressurisation
Many parameters can be automated
Carbonates foaming products.

Technical Specification :

Carbonator Filler

Laboratory Instrument India

