



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Carbonator Filler

**Product Code :**  
LBNY-0005-10700031



**Description :**

Accurate control of filling and carbonating variables  
Throughputs up to 60l/h achievable (2l bottles)  
Wide range of cans or bottles can be used from 0.15l to in excess of 2l  
Evaluation service for customer cans and bottles  
Option of filling and crown seal capping and/or screw top capping at a single station  
15l process vessel, positioned at a convenient height  
Batch & continuous processing  
CO<sub>2</sub> content in excess of 10g/l can be realised  
Enhanced colour touch screen control  
Carbonates at any temperature between 4°C and 20°C  
Premix and postmix capability  
Fills non-carbonated drinks  
In-bottle pressure measurement enables accurate and reproducible depressurisation  
Many parameters can be automated  
Carbonates foaming products.

**Technical Specification :**

Carbonator Filler

---

# Laboratory Instrument India

