



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Laboratory Scale Rapid Extractor

**Product Code :**  
LBNY-0005-10700032



**Description :**

High pressure, room temperature extraction, excellent for sensitive materials.  
Easy to use: microprocessor control and real time display of sequencing and process parameters.  
Fast, efficient solid-liquid extraction of active constituents from plants, herbs, fruits and other similar materials.  
Unaltered natural properties and characteristics of the active constituents, no discolouration or oxidisation  
Automatic alarm indication of process anomalies  
Highly flexible extraction parameters.  
Automatic pressure control  
Automatic solvent loading and priming.

**Technical Specification :**

Laboratory Scale Rapid Extractor



**Laboratory Instrument India**