



Email : sales@laboratoryinstrumentindia.com

Product Name :
Laboratory Scale Rapid Extractor

Product Code :
LBNY-0005-10700032



Description :

High pressure, room temperature extraction, excellent for sensitive materials.
Easy to use: microprocessor control and real time display of sequencing and process parameters.
Fast, efficient solid-liquid extraction of active constituents from plants, herbs, fruits and other similar materials.
Unaltered natural properties and characteristics of the active constituents, no discolouration or oxidation
Automatic alarm indication of process anomalies
Highly flexible extraction parameters.
Automatic pressure control
Automatic solvent loading and priming.

Technical Specification :

Laboratory Scale Rapid Extractor



Laboratory Instrument India