



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Benchtop Rapid Extractor

**Product Code :**  
LBNY-0005-10700036



**Description :**

Unaltered natural properties and characteristics of the active constituents, no discolouration or oxidation  
High pressure, room temperature extraction, excellent for sensitive materials  
Easy to use: microprocessor control and real time display of sequencing and process parameters  
Fast, efficient solid-liquid extraction of active constituents from plants, herbs, fruits and other similar materials  
Automatic pressure control  
Highly flexible extraction parameters  
Automatic alarm indication of process anomalies.

**Technical Specification :**

Benchtop Rapid Extractor



**Laboratory Instrument India**