



Email : sales@laboratoryinstrumentindia.com

Product Name :
Benchtop Rapid Extractor

Product Code :
LBNY-0005-10700036



Description :

Unaltered natural properties and characteristics of the active constituents, no discolouration or oxidation
High pressure, room temperature extraction, excellent for sensitive materials
Easy to use: microprocessor control and real time display of sequencing and process parameters
Fast, efficient solid-liquid extraction of active constituents from plants, herbs, fruits and other similar materials
Automatic pressure control
Highly flexible extraction parameters
Automatic alarm indication of process anomalies.

Technical Specification :

Benchtop Rapid Extractor



Laboratory Instrument India