



Email : sales@laboratoryinstrumentindia.com

Product Name :
Neutraliser Washer Bleacher

Product Code :
LBNY-0005-10700037



Description :

Agitation by variable speed flat blade impeller (variable to 1800rpm)
Working vacuum of 75mm Hg (ABS) with liquid ring vacuum pump
Simple Controls
Surfaces in contact with the process fluid are stainless steel
Compact and self contained
Cooling coil
Pressure leaf filter designed
25 litre vessel designed
Electrical heating element
Extraction pump.

Technical Specification :

Neutraliser Washer Bleacher



Laboratory Instrument India

