



Email : sales@laboratoryinstrumentindia.com

Product Name :
Sterile Filling System

Product Code :
LBNY-0005-10700038



Description :

Capable of filling containers and producing sterile samples
Simple operation
Large working area
Open front for ease of use
Filling environment
Wide range of containers and sealing devices can be catered for
Cleaning, Sterilisation and Filling modes
Stainless Steel Construction
Low Oxygen Filling facility as standard
Can form part of a complete aseptic UHT processing line.

Technical Specification :

Sterile Filling System



Laboratory Instrument India

