



Email : sales@laboratoryinstrumentindia.com

Product Name :
Deodorising Unit

Product Code :
LBNY-0005-10700039



Description :

10 litre capacity stainless steel jacketed vat
Steam ejector
Condenser
Stainless steel pressure vessel
Electrical heating of the oil
Comprehensive control panel and instructions
Liquid ring vacuum pump
Internal cooling coil
Polishing filter
Extraction pump
Vacuum to 1 Torr
Internal steam sparge.
Temperatures to 250°C.

Technical Specification :

Deodorising Unit

Laboratory Instrument India

