



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Multi Purpose Processing Vessel

**Product Code :**  
LBNY-0005-10700043



**Description :**

The unit is fully self-contained with electrical heating, cooling water and chilled water circulating systems and an efficient emulsifier/mixer, which is easily raised or lowered into the vat.

A comprehensive control console, also housed in the cabinet, incorporates a chart recorder which continuously records the temperature of the vat contents. Also included on the console is a digital readout for salient temperatures and heating controls.

The consists of a jacketed vat, mounted within a mobile cabinet containing all the servicing utilities and necessary instrumentation.

Hygienic interior of vessel enlarge

Hygienic interior of vessel

Schematic diagram of the Schematic diagram of the enlarge

The vat is constructed from high quality stainless steel, the interior having a mirror finish.

**Technical Specification :**

Multi Purpose Processing Vessel



**Laboratory Instrument India**

