



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Oil Extraction Screw Press

**Product Code :**  
LBNY-0005-10700045



**Description :**

The press is driven at the discharge end through a single reduction helical gearbox. Seed or ground meal is gravity fed to the press from an integral feed hopper fitted with a manually operated shut-off slide. The power source is a 4kW, 4 pole, 3 phase electric motor. The thickness of the press cake may be varied.

**Technical Specification :**

**Technical Details:**

Complete with feed hopper, manual feed chute, oil discharge chute and suitable base plate. Includes: one 4kW, 4 pole motor. Screw press, having barrel constructed from separate cast rings spaced apart by shims. 60mm diameter bore, 155mm drainage length with single piece wormshaft driven at the discharge end through a single reduction helical gearbox.



**Laboratory Instrument India**