



Email : sales@laboratoryinstrumentindia.com

Product Name :
Rising Film Evaporator

Product Code :
LBNY-0005-10700048



Description :

The system is supported within a stainless steel framework and requires only connection to suitable steam, electrical and water services for operation.

The component parts of the industrial process have been scaled down, feed tank, vapour / liquid separator, condenser, concentrate tank, condensate tank, CIP pump and vacuum pump, all manufactured from hygienic materials such as stainless steel, PTFE and glass. The single tube-in-jacket evaporator replicates industrial units which often contain scores of much longer tubes within the steam space.

Only 0.65m² of floor space is required to install the evaporator and it has been designed to fit easily into a room of standard height.

Technical Specification :

Rising Film Evaporator



Laboratory Instrument India