




Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Plate And Frame Filter Press

**Product Code :**  
LBNY-0005-10700051



**Description :**

A centrifugal pump is mounted integrally in the framework and is connected to the filter inlet via a throttle valve and pressure gauge.

A fully self-contained plate and frame type filter press constructed entirely in stainless steel.

Filter surface area is 0.22m<sup>2</sup> and solids holding capacity is 1.5 litres.

A back pressure valve is located at the filter outlet.

The filter pack is made up of four plates and four frames, which will accept a dressing of seven filter sheets 200mm x 200mm.

The filter pack and sheets are compressed manually using a hardened steel screw thread.

**Technical Specification :**

Plate And Frame Filter Press



**Laboratory Instrument India**