



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Autoclave

**Product Code :**  
LBNY-0005-10700052



**Description :**

Die-cast in special high tensile aluminium alloy with an inner aluminium liner to facilitate cleaning. The heavy duty lid is secured by means of butterfly swivel clamps, two of which act as safety retaining bolts.

A compact pressure steriliser of high quality for use at 0.75 or 1.05kg/cm<sup>3</sup>, complete with built-in heating elements, control switches and fitted with resettable microtherm cut-out.

Supplied complete with adjustable safety valve, pressure gauge, air cock, drain off syphon and inner liner. Also supplied, 10 robust autoclave milk sample bottles and screw caps, 250ml capacity, and 55 metres of autoclave tape.

**Technical Specification :**

Autoclave



**Laboratory Instrument India**