



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Can Bus In The Comfort System

**Product Code :**  
LBNY-0005-10800024



#### **Description :**

The stand allows the performing of a diagnosis of the system through the OBD2 connector that is connected to the main comfort controller. Complete with user manual that will allow teachers to develop their own experiments. The electric input and output signals are transmitted through banana plugs for quick measurement of their parameters. There is the possibility to change the lock and alarm mode through a new controller coding. This demonstration panel presents the functioning of the comfort system based on data transfer performed with the aid of CAN BUS.

The Can Bus In The Comfort Systems consists of:  
Electric motors used to wind down and up the glasses  
Electric door locks  
Heated side mirrors equipped with position adjustment mechanism  
A set of switches for the control of the actuators  
Alarm  
The comfort system controller  
Simulators of hood and tailboard locks  
Vehicle interior lighting.

#### **Technical Specification :**

Can Bus In The Comfort System



**Laboratory Instrument India**