



Email : sales@laboratoryinstrumentindia.com

Product Name :
Vehicle Sensor System

Product Code :
LBNY-0005-10800025



Description :

Vehicle Sensor System demonstration panel enables learning in terms of combination, test and evaluation of the parameters of the Vehicle Sensor System sub-assemblies. It can be configured in different ways.

The Vehicle Sensor System trainer includes the following systems:

- Engine and air temperature sensor
- Lambda sensor
- Rotation speed active sensor
- Mass and volume air flow meters check system
- Knock sensor
- Set of the main vehicle system sensors
- Differential pressure sensor
- Speedometer sensor
- Rotation direction sensor
- Fuel level sensor
- Oil pressure sensor
- Acceleration sensor.

Technical Specification :

Vehicle Sensor System



Laboratory Instrument India