




Email : sales@laboratoryinstrumentindia.com

Product Name :
Autotronic Training System

Product Code :
LBNY-0005-10800031



Description :

The whole system has been modularized to provide electronic and automotive courses with step-by-step experiments on technological educational training.

Interactive Computerized Automotive System

The simulate the operation of fuel injection system, ignition system, and the control of exhaust gas, etc.

Experiments include the sensors characteristic and operation of various sensors and actuators monitored by microprocessor on the main unit.

All theoretical, experimental and practical learning procedures are supported by a personal computer assisted supervision and dedicated software.

Technical Specification :

Autotronic Training System



Laboratory Instrument India