



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Hydraulics Bench Chronometer Digital Flow Meter

**Product Code :**  
LBNY-0002-710030



**Description :**

Hydraulics Bench Chronometer Digital Flow Meter

**Technical Specification :**

Hydraulics Bench, Chronometer, Digital Flow Meter

The Hydraulics Bench, allows the study of fluid behaviour, hydraulic theory and the properties of fluid mechanics. Base unit to provide water supply and volumetric measurement to many other experimental accessories.

Molded fiberglass bench-top features an open channel and a measuring tank with space for mounting accessories. Circulating pump.

Fiberglass main storage tank with drain valve. Tank is stepped for measurement of both high and low flow rates.

Level indicators for water levels in the open channel and measuring tank.

Bull's eye level for correct level reading. Variable area flow meter (optional). Mounted on locking wheels..

Hydraulics Bench Chronometer Digital Flow Meter, Hydraulics Bench Chronometer Digital Flow Meter Equipments, Hydraulics Bench Chronometer Digital Flow Meter Tools, Hydraulics Bench Chronometer Digital Flow Meter Tool Kits, Hydraulics Bench Chronometer Digital Flow Meter Manufacturers, Hydraulics Bench Chronometer Digital Flow Meter Suppliers from India, China, Kenya

---

# Laboratory Instrument India

