



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Instructional Kit, Assembly with Taper Keys

**Product Code :**  
LBNY-0005-10900016



**Description :**

This instructional kit provides the necessary material for students to learn systematically how a shaft is joined to a hub or a coupler in a professional way using taper keys. The material is clearly laid out on a plastic tray. We recommend that the exercises are carried out in a workshop, as all preparations such as filing, drilling, countersinking and tapping must be carried out in a correct and proper manner by the student. The well-structured instructional material outlines all the necessary technical information and provides a step-by-step guide through the exercises.

**Technical Specification :**

Set of material for workshop exercises in taper key joints  
Key driver, key extractor, driver sleeve  
1 Set of semi-finished items for the production of taper keys to DIN 6886 and DIN 6887  
Shaft with slots  
2 Hubs with slots  
Coupler  
Multiple trays stackable  
All parts clearly laid out on a tray.

---

# Laboratory Instrument India

