



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Drilling Kit

**Product Code :**  
LBNY-0005-10900017



**Description :**

The kit is used primarily for viewing and information purposes. The cutting geometry has been deliberately changed on some of the drills so that the influence of the cutting angle and clearance angle can be demonstrated. Exercises or experiments can only be carried out to a limited extent. It includes 19 different drills, including special types, e.g. a centre drill or a taper pin hole drill. The well-structured instructional material enhances the informational value of the collection. In addition a taper shaft and the associated key, as well as a sliding gauge are included. The kit is clearly laid out on a plastic tray.

**Technical Specification :**

Drilling tools demonstration collection  
All parts clearly laid out on a plastic tray  
Multiple trays stackable  
Content: 19 different drills, 7 workpieces with specimen cuts, 1 taper shaft MK1 with key and a twist drills grinding gauge for checking point angle.



**Laboratory Instrument India**

