



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Milling Kit

**Product Code :**  
LBNY-0005-10900020



**Description :**

The kit is used primarily for viewing and information purposes. No provision is made for conducting exercises or experiments. The well-structured instructional material enhances the informational value of the kit. It includes 12 different milling cutters. A shell end mill arbor and 2 cutter arbor rings are also included. A small steel plate with examples of machining is provided to help identify the possible applications for the different types of milling cutter. The kit is clearly laid out on a plastic tray.

**Technical Specification :**

Milling cutter demonstration kit

Multiple trays stackable

Content: 12 different milling cutters, 1 shell end mill arbor with 2 cutter arbor rings and 1 steel plate with examples of machining.

All parts clearly laid out on a plastic tray.



**Laboratory Instrument India**