



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Countersinking Kit

**Product Code :**  
LBNY-0005-10900021



#### **Description :**

The kit is used primarily for viewing and information purposes. No provision is made for conducting exercises or experiments. The largest tool diameter is 16,75mm. The kit is clearly laid out on a plastic tray. The well-structured instructional material enhances the informational value of the collection. It also includes interchangeable guiding pins and a machining specimen with initial countersinking. It includes 12 different countersinks including special types, e.g. a backward countersink tool.

#### **Technical Specification :**

Countersinking tools demonstration collection

Multiple trays stackable

All parts clearly laid out on a plastic tray

Content: 12 different countersink tools, 3 interchangeable guiding pins, 1 holder for a backward countersink tool and a specimen with initial countersinking.

#### **Technical Data:**

1 Core drill: D=16,75mm

1 Backward countersink tool: D=15mm

1 Counterbore: D=15mm

5 Conical countersink tools (M8/90°, C20/60°, A20/60°, C15/90°, C16, 5/90°)

4 Flat countersink tools (3xM8, D=15mm).



**Laboratory Instrument India**