



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Turning Kit

**Product Code :**  
LBNY-0005-10900022



**Description :**

The kit is used primarily for viewing and information purposes. The kit is clearly laid out on a plastic tray. The well-structured instructional material enhances the informational value of the kit. No provision is made for conducting exercises or experiments. Four reversible carbide tips and a turned part are also included. The turned part features examples of machining with the different tools. It includes 13 different lathe tools, including a reversible carbide tip holder.

**Technical Specification :**

Lathe tools demonstration kit  
Multiple trays stackable  
Content: 13 different lathe tools, 4 reversible carbide tips, 1 turned part  
Turned part: aluminium shaft 262mm long, diameter 70mm.  
All parts clearly laid out on a plastic tray.

**Technical Data:**

1 Reversible carbide tip holder  
12 Lathe tools with soldered-in carbide tips (including Inside tool, parting-off tool, side tool, face tool)  
4 Reversible carbide tips of different materials.

---

# Laboratory Instrument India

