



Email : sales@laboratoryinstrumentindia.com

Product Name :
Turning Kit

Product Code :
LBNY-0005-10900022



Description :

The kit is used primarily for viewing and information purposes. The kit is clearly laid out on a plastic tray. The well-structured instructional material enhances the informational value of the kit. No provision is made for conducting exercises or experiments. Four reversible carbide tips and a turned part are also included. The turned part features examples of machining with the different tools. It includes 13 different lathe tools, including a reversible carbide tip holder.

Technical Specification :

Lathe tools demonstration kit
Multiple trays stackable
Content: 13 different lathe tools, 4 reversible carbide tips, 1 turned part
Turned part: aluminium shaft 262mm long, diameter 70mm.
All parts clearly laid out on a plastic tray.

Technical Data:

1 Reversible carbide tip holder
12 Lathe tools with soldered-in carbide tips (including Inside tool, parting-off tool, side tool, face tool)
4 Reversible carbide tips of different materials.

Laboratory Instrument India

