



Email : sales@laboratoryinstrumentindia.com

Product Name :
Reaming Kit

Product Code :
LBNY-0005-10900023



Description :

The kit is used primarily for demonstration and information purposes. It includes 10 different reaming tools and a limit plug gauge. It allows only limited scope for conducting exercises or experiments. The well-structured instructional material enhances the informational value of the kit. The kit includes a sample plate with reamed bores with which a fit can be checked. The kit is clearly laid out on a plastic tray.

Technical Specification :

Reaming tools demonstration collection

Multiple trays stackable

All parts clearly laid out on a plastic tray

Content: 10 different reamers, 1 limit plug gauge and a sample plate with reamed bores.

Technical Data:

2 Hand reamers, 1 of which is adjustable, D=10mm

2 Machine reamers with Morse taper: D=5mm (Taper 1:50) and MK1

Limit plug gauge: D=10mm, fit H7

1 Stub reamer: D=10mm

7 Chucking reamers: D=10mm.

Laboratory Instrument India

