



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Reaming Kit

**Product Code :**  
LBNY-0005-10900023



#### **Description :**

The kit is used primarily for demonstration and information purposes. It includes 10 different reaming tools and a limit plug gauge. It allows only limited scope for conducting exercises or experiments. The well-structured instructional material enhances the informational value of the kit. The kit includes a sample plate with reamed bores with which a fit can be checked. The kit is clearly laid out on a plastic tray.

#### **Technical Specification :**

Reaming tools demonstration collection  
Multiple trays stackable  
All parts clearly laid out on a plastic tray  
Content: 10 different reamers, 1 limit plug gauge and a sample plate with reamed bores.

#### **Technical Data:**

2 Hand reamers, 1 of which is adjustable, D=10mm  
2 Machine reamers with Morse taper: D=5mm (Taper 1:50) and MK1  
Limit plug gauge: D=10mm, fit H7  
1 Stub reamer: D=10mm  
7 Chucking reamers: D=10mm.

---

# Laboratory Instrument India

