



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Dimensional Metrology I Training Kit One

**Product Code :**  
LBNY-0005-10900030



**Description :**

The dimensional metrology practice kit is designed for practical training in the metalworking trades. They were manufactured with the accuracy of CNC parts. All parts are clearly laid out and well protected on a plastic tray. The trays are stackable, providing for space-saving storage. The kit includes ten spacer plates, all of which differ to a minor degree in length and each of which is individually marked. The kit offers the advantage of being ready for use immediately; it contains everything required for the exercises. The test pieces used are flat stainless steel spacer plates.

**Technical Specification :**

All parts clearly laid out and well protected on a plastic tray  
Measurement exercises on a stainless steel spacer plate  
Instructional kit complete with test pieces and measuring aids  
10 test pieces, each of different dimensions  
Practice kit for dimensional metrology in the metalworking trades

**Technical Data:**

Test pieces: 13 lengths measurable  
Steel ruler: 300mm  
Vernier caliper gauge: 200mm  
Set square: 90, l x w: 100x70mm.



**Laboratory Instrument India**