



Email : sales@laboratoryinstrumentindia.com

Product Name :
Cutaway Model Planetary Gear

Product Code :
LBNY-0005-10900032



Description :

The range consists of a series of different cutaway models of drive system components. The mechanical functionality of all components is retained in full. Each of the cutaway models is securely mounted on a base plate to which handles are attached for carrying purposes. All mechanisms are hand-driven. All components are industrially relevant and have been prepared such that individual active elements are clearly visible to the user.

Technical Specification :

Industrially-relevant part, fully functional cutaway model
Solid metal base plate, with handles
Hand-operated tabletop cutaway model for the demonstration of the functionality of a planetary gear
Dimensions and Weight:
l x w x h : 350 x 300 x 160 mm
Weight: 5kg.



Laboratory Instrument India