



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Cutaway Model Hermetic Refrigerant Compressor

**Product Code :**  
LBNY-0005-10900034



**Description :**

Hermetic refrigerant compressors have welded sheet steel housing. The drive motor and refrigerant compressor are fitted in this capsule. The motor compressor is suspended by springs in the capsular housing which minimises the transmission of vibration to the capsule body, thereby reducing noise emission. Unlike open compressors, the refrigerant compressor is directly connected to the drive motor and therefore does not require any floating seal.

**Technical Specification :**

**Technical Data:**

1-Cylinder compressor according to CECOMAF

Power consumption: 321W

Refrigerating capacity at 5/55°C: 592W

Dimensions and Weight:

LxWxH: 300x300x240mm

Weight: 12kg.



**Laboratory Instrument India**

