



Email : sales@laboratoryinstrumentindia.com

Product Name :
Dimensional Metrology I Training Kit Five

Product Code :
LBNY-0005-10900035



www.laboratoryinstrumentindia.com

Description :

The kit offers the advantage of being ready for use immediately; it contains everything required for the exercises. This dimensional metrology practice kit is designed for practical training in the metalworking trades. All parts are clearly laid out and well protected on a plastic tray. The kit includes ten shafts, all of which differ to a minor degree in dimensions and each of which is individually marked. They were manufactured with the accuracy of CNC parts. The trays are stackable, providing for space-saving storage. The test pieces used are stainless steel shafts with various protrusions.

Technical Specification :

Measurement exercises on a stainless steel shaft
Instructional kit complete with test pieces and measuring aids
Practice kit for dimensional metrology in the metalworking trades
10 Test pieces, each of different dimensions
Detailed instructional material
All parts clearly laid out and well protected on a plastic tray.



Laboratory Instrument India

