




Email : sales@laboratoryinstrumentindia.com

Product Name :
Cutaway Model Spur Gear

Product Code :
LBNY-0005-10900036



Description :

The range consists of a series of different cutaway models of drive system components. The mechanical functionality of all components is retained in full. Each of the cutaway models is securely mounted on a base plate to which handles are attached for carrying purposes. All components are industrially relevant and have been prepared such that individual active elements are clearly visible to the user. All mechanisms are hand-driven.

Technical Specification :

Industrially-relevant part, fully functional cutaway model.
Solid metal base plate, with handles.
Hand-operated tabletop cutaway model for the demonstration of the functionality of a spur gear.



Laboratory Instrument India