



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Drilling Jig For Flat Part

**Product Code :**  
LBNY-0005-10900037



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

The surfaces of the steel parts have a gunmetal finish to prevent corrosion. Comprehensive and well-structured instructional material supports wide-ranging use of the model in technical training. The drilling jig model forms part of the extensive course providing an introduction to technical drawing. The approach involves teaching the systematic and transferable learning matter based on an operational mechanism. All parts are precision-manufactured on CNC machines.

**Technical Specification :**

All drilling jig components made from steel, precision-manufactured, with gunmetal surface finish. Instructional material incorporates action-oriented and interdisciplinary forms of teaching.  
Part of the technical drawing course.  
Practical drilling jig for machining of a planar part.  
Dimensions and Weight:  
l x w x h: 70 x 60 x 8 mm (plane part)  
l x w x h: 90 x 90 x 40 mm (drilling jig)  
Weight: 2 kg



**Laboratory Instrument India**

