



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Dimensional Metrology I Training Kit Two

**Product Code :**  
LBNY-0005-10900038



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

#### **Description :**

This dimensional metrology practice kit is designed for practical training in the metalworking trades. The kit includes ten bearing plates, all of which differ to a minor degree in dimensions and each of which is individually marked. The test pieces used are stainless steel bearing plates. They were manufactured with the accuracy of CNC parts. The kit offers the advantage of being ready for use immediately; it contains everything required for the exercises. The trays are stackable, providing for space-saving storage. All parts are clearly laid out and well protected on a plastic tray.

#### **Technical Specification :**

All parts clearly laid out and well protected on a plastic tray  
Training kit for dimensional metrology in the metalworking trades  
Instructional kit complete with test pieces and measuring aids  
Measurement exercises on a stainless steel bearing plate  
10 Test pieces, each of different dimensions

#### **Technical Data:**

Vernier caliper gauge: 0...200mm  
External micrometer: 0...25mm, resolution: 0,01mm  
Depth micrometer: 0...25mm  
Pocket caliper gauge: 0...150mm  
Depth caliper gauge: 0...150mm  
Test pieces  
9 Lengths, 4 depths, 4 diameters measurable

---

Dimensions and Weight:

l x w x h: 500x350x110mm (tray)

Weight: 7kg.



**Laboratory Instrument India**