



Email : sales@laboratoryinstrumentindia.com

Product Name :
Dimensional Metrology II Training Kit Two

Product Code :
LBNY-0005-10900039



www.laboratoryinstrumentindia.com

Description :

This dimensional metrology practice kit is designed for practical training in the metalworking trades. They were manufactured with the accuracy of CNC parts. The kit includes eight tapers, all of which differ to a minor degree in dimensions and each of which is individually marked. The kit offers the advantage of being ready for use immediately; it contains everything required for the exercises. The test pieces used are stainless steel internal and external tapers.

Technical Specification :

Instructional kit complete with test pieces and measuring aids
8 Stainless steel test pieces, each of different dimensions
Practice kit for dimensional metrology in the metalworking trades
Measurement exercises on an external and internal taper
All parts clearly laid out and well protected on a plastic tray

Technical Data:

Depth caliper gauge: 0...150mm
Test pieces: 3 dimensions measurable
Taper plug gauge MK 3
Taper ring gauge MK 3
Dimensions and Weight:
LxWxH: 500x350x110mm (tray)
Weight: approx. 7kg.



Laboratory Instrument India