



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Dimensional Metrology I Training Kit Six

**Product Code :**  
LBNY-0005-10900040



#### **Description :**

This dimensional metrology practice kit is designed for practical training in the metalworking trades. The test pieces used are stainless steel shafts with various protrusions including one square section. The kit includes ten shafts, all of which differ to a minor degree in dimensions and each of which is individually marked. The kit offers the advantage of being ready for use immediately; it contains everything required for the exercises. All parts are clearly laid out and well protected on a plastic tray. They were manufactured with the accuracy of CNC parts. The trays are stackable, providing for space-saving storage.

#### **Technical Specification :**

Instructional kit complete with test pieces and measuring aids  
All parts clearly laid out and well protected on a plastic tray  
10 Test pieces, each of different dimensions  
Practice kit for dimensional metrology in the metalworking trades  
Measurement exercises on a stainless steel shaft

#### **Technical Data:**

Radius gauges: 1...7mm concave and convex  
Vernier caliper gauge: 0...200mm  
Depth caliper gauge: 0...150mm  
Test pieces  
9 Lengths and angles, 7 diameters and radii measurable.



**Laboratory Instrument India**