



Email : sales@laboratoryinstrumentindia.com

Product Name :
Insulated Torque Wrench

Product Code :
LBNY-0002-660007



Description :

Insulated Torque Wrench

Technical Specification :

Suitable for control of inspection and test equipment.
Sustainable and environmentally sound.
Fine scale graduation with vernier scale.
Locking against unintentional adjustment of the setting value.
(Torque Wrench With Protective Insulation 1000V)

Insulated Torque Wrench, Insulated Torque Wrench Equipments, Insulated Torque Wrench Tools, Insulated Torque Wrench Tool Kits, Insulated Torque Wrench Manufacturers, Insulated Torque Wrench Suppliers from India, China, Kenya



Laboratory Instrument India