



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Insulated Torque Wrench

**Product Code :**  
LBNY-0002-660007



**Description :**

Insulated Torque Wrench

**Technical Specification :**

Suitable for control of inspection and test equipment.  
Sustainable and environmentally sound.  
Fine scale graduation with vernier scale.  
Locking against unintentional adjustment of the setting value.  
(Torque Wrench With Protective Insulation 1000V)

Insulated Torque Wrench, Insulated Torque Wrench Equipments, Insulated Torque Wrench Tools, Insulated Torque Wrench Tool Kits, Insulated Torque Wrench Manufacturers, Insulated Torque Wrench Suppliers from India, China, Kenya



**Laboratory Instrument India**