



Email : sales@laboratoryinstrumentindia.com

Product Name :
Dimensional Metrology I Training Kit Four

Product Code :
LBNY-0005-10900041



www.laboratoryinstrumentindia.com

Description :

This dimensional metrology practice kit is designed for practical training in the metalworking trades. The test pieces used are stainless steel angle plates that have been manufactured with the accuracy of CNC parts. The kit includes ten angle pieces, all of which differ to a minor degree in dimensions and each of which is individually marked. The kit offers the advantage of being ready for use immediately; it contains everything required for the exercises. The trays are stackable, providing for space-saving storage. All parts are clearly laid out and well protected on a plastic tray.

Technical Specification :

Measurement exercises on a stainless steel angle piece
Instructional kit complete with test pieces and measuring aids
All parts clearly laid out and well protected on a plastic tray
Practice kit for dimensional metrology in the metalworking trades
10 Test pieces, each of different dimensions.



Laboratory Instrument India

