



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Dimensional Metrology II Training Kit One

**Product Code :**  
LBNY-0005-10900042



#### **Description :**

This dimensional metrology practice kit is designed for practical training in the metalworking trades. The test pieces used are stainless steel slotted shafts. They were manufactured with the accuracy of CNC parts. The kit offers the advantage of being ready for use immediately; it contains everything required for the exercises. All parts are clearly laid out and well protected on a plastic tray. The trays are stackable, providing for space-saving storage. The kit includes six shafts, all of which differ to a minor degree in dimensions and each of which is individually marked.

#### **Technical Specification :**

Instructional kit complete with test pieces and measuring aids  
All parts clearly laid out and well protected on a plastic tray  
6 Test pieces, each of different dimensions  
Practice kit for dimensional metrology in the metalworking trades  
Measurement exercises on a stainless steel shaft

#### **Technical Data:**

Test pieces: 4 dimensions measurable  
Dial gauge: 0...1mm, graduations: 0,001mm  
Key-way vernier caliper: D=5...80mm  
V-Block: indent 90  
Pocket caliper gauge: 0...150mm  
Depth micrometer: 0...75mm  
Dimensions and Weight:

---

LxWxH: 500x350x110mm (tray)  
Weight: 12kg.



**Laboratory Instrument India**