



Email : sales@laboratoryinstrumentindia.com

Product Name :
Cylindrical Work Samples With Slanted Cut-Outs

Product Code :
LBNY-0005-10900043



www.laboratoryinstrumentindia.com

Description :

The set contains 18 aluminium drawing models. The cylindrical models have cut-outs parallel to, and at an angle to, the main axis. The models are precision-manufactured and are suitable for measurement exercises as well as drawing exercises. The geometric models support the learning process by providing an introduction to technical drawings: from the solid model to the more abstract representation of the three views in a technical drawing. Two students can easily work with one model set. For optimal learning, it is recommended that students work independently from each other on these exercises.

Technical Specification :

Finely graduated degrees of difficulty within the model series
All parts clearly laid out on a practical tray
Model set providing an introduction to technical drawing
Set of 18 cylindrical models, made from aluminium, with anodised finish
Models with slanted cut-outs
Multiple trays stackable
Precision manufacture of the models enables measurement exercises
Lesson-supporting instructional material

Technical Data:

Dimensions and Weight of Cylindrical Work Samples with Slanted Cut-Outs
l x w x h: 600 x 400 x 140 mm (tray)
Dimensions of each model: d x h: 50x80mm
Weight: 9 kg.



Laboratory Instrument India