



Email : sales@laboratoryinstrumentindia.com

Product Name :
Drilling Jig For An Annular Disc

Product Code :
LBNY-0005-10900045



Description :

All parts, made from steel, are precision-manufactured on CNC machines. The drilling jig model forms part of the extensive course providing an introduction to technical drawing. The approach involves teaching the systematic and transferable learning matter based on an operational mechanism. Comprehensive and well-structured instructional material supports wide-ranging use of the model in technical training. The surfaces of the steel parts have a gunmetal finish to prevent corrosion.

Technical Specification :

All drilling jig components made from steel, precision-manufactured, with gunmetal surface finish
Instructional material incorporates action-oriented and interdisciplinary forms of teaching
Part of the technical drawing course
Practical drilling jig for machining of a bearing cover
Dimensions and Weight:
D x h: 90 x 10 mm (ring)
D x h: 120 x 70 mm (drilling jig)
Weight: 3 kg.



Laboratory Instrument India

