



Email : sales@laboratoryinstrumentindia.com

Product Name :
Lever Shears

Product Code :
LBNY-0005-10900046



www.laboratoryinstrumentindia.com

Description :

The lever shears model can be considered to be a project enabling interdisciplinary, action-oriented teaching. As well as the primary technical drawing focus of teaching, other topics that can be covered include machine elements, assembly processes, and manufacturing technology. The lever shears model forms part of the extensive course providing an introduction to technical drawing. Comprehensive and well-structured instructional material supports wide-ranging use of the model in technical training. The approach involves teaching the systematic and transferable learning matter based on an operational mechanism.

Technical Specification :

Instructional material incorporates action-oriented and interdisciplinary forms of teaching
All lever shears components made from steel, precision-manufactured, with gunmetal surface finish
Complex function model of a lever shears
Part of the technical drawing course.
Dimensions and Weight of Lever Shears
l x w x h: 120x120x120mm (excluding lever)
Weight: 3 kg.



Laboratory Instrument India

