



Email : sales@laboratoryinstrumentindia.com

Product Name :
Stirred Tank Reactors In Series

Product Code :
LBNY-0005-1100003



Description :

A dead-time residence coil can also be attached to the exit of the last reactor in the series.
A self-contained bench mounted small scale unit fitted with three continuous stirred reactors in series which are fed from two 5 litre tanks. Each reactor is fitted with a conductivity probe.
There are two independent, variable speed feed pumps.

Technical Specification :

Stirred Tank Reactors In Series



Laboratory Instrument India