



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Stirred Tank Reactors In Series

**Product Code :**  
LBNY-0005-1100003



**Description :**

A dead-time residence coil can also be attached to the exit of the last reactor in the series.  
A self-contained bench mounted small scale unit fitted with three continuous stirred reactors in series which are fed from two 5 litre tanks. Each reactor is fitted with a conductivity probe.  
There are two independent, variable speed feed pumps.

**Technical Specification :**

Stirred Tank Reactors In Series



**Laboratory Instrument India**