



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Computer Controlled Tray Dryer

**Product Code :**  
LBNY-0005-1100006



**Description :**

A compact desktop tray dryer suitable for laboratory use.  
Computer control of temperature and air flow rate, with sophisticated data logging and analysis software.  
Integrated electronic weight measurement to determine drying rate.  
Electronic measurement of air flow.  
Stainless Steel Construction.  
Electronic measurement of temperature and humidity before and after the drying trays.  
Capacity: up to 2.1 Kg of wet material.  
Flow rates, 0.4 - 3.0m/s over trays.  
Temperatures up to 80°C at 0.4 m/s less at higher flows.

**Technical Specification :**

Computer Controlled Tray Dryer



**Laboratory Instrument India**